

# Vibration Sensitivity of a Laboratory Bench Microscope

Hal Amick

*Colin Gordon & Associates, USA*

Matthew Stead

*Bassett Acoustics, Australia*

Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the  
Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Overview

- Definition of Problem
- Tests of Microscope Base Sensitivity
- Floor Vibrations
- Microscope Support
  - Common Types
    - Laboratory Benches
    - Optical
  - Vibrations at Base of Microscope
- Vibration Criteria as a Function of Magnification and Support


Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Definition of Problem

- Vibration degrades image observed in microscope
  - Low frequency: Jiggle
  - Higher frequency: Blurring
- Not much information regarding vibration sensitivity of optical microscopes
- How does sensitivity vary with magnification?
- How the bench contribute?
- How much does an isolation table help?

Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Typical Walker-Generated Floor Vibration (Vertical)



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Tests of Microscope Base Sensitivity

- Benchtop represents interface between the microscope and “the rest of the world”
- Document the sensitivity at the microscope’s base
- Document the modification to floor motion provided by the support bench
- “Back out” the limits of floor vibration

Hal Amick & Michael Stead, “Vibration Sensitivity of a Laboratory Bench Microscope,” Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Microscopes Used in Study

- Nikon Eclipse E400
  - 40x, 100x, 400x
- Nikon Optiphot
  - 1000x



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Basic Test Configuration

- Microscope and shaker on “springy” table
- Graduated Calibration Standard for gauging motion
- Tested in three directions:
  - *Vertical* – Microscope placed midspan, shaker placed behind it, vertical force
  - *Fore-and-aft* – Microscope placed at one end on centerline, shaker placed on table behind it with force axis through microscope
  - *Side-to-side* – Microscope placed at middle of edge on long side, shaker on table behind it above long’l axis of table

# Test a Given Direction's Sensitivity


1. Select a frequency
2. Set at lowest amplitude, turn on shaker
3. Watch through eyepiece
4. Increase amplitude until motion first visible
5. Record frequency and amplitude
6. Change frequency, repeat

Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico



# Lower Bounds – 1


*(Omnidirectional)*



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Lower Bounds – 2

*(Omnidirectional)*



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Microscope Support


- Primary excitation is floor motion due to people walking
  - Laboratory – Predominantly vertical
  - Raised access floor – Predominantly horizontal
- On laboratory bench
  - Bench amplifies vertical motion; “stiff” horizontally
- On cleanroom bench
  - Bench amplifies both horizontal and vertical motion
- On pneumatically isolated bench or benchtop support, such as those from Newport Corporation
  - Bench amplifies at very low frequencies and attenuates at higher frequencies

# Field Tests

- Performed at UCSF Genetech Hall, biotechnology research lab under construction
  - Built-in laboratory bench
  - Newport microscope isolation bench
- Measurement of walker-generated vertical vibrations on floor and benchtops
- Measurement of transmissibility of lab bench and Newport bench


Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Compare Transmissibility (Vertical Only)




Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Time Histories Floor & Newport Table




Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Time Histories Floor & Newport Table



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Walker-Generated Spectra Floor & Newport Table



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico




# Vibration Criteria

1/3 O.B. Amplitude, $\mu\text{m/s}$ ( $\mu\text{in/s}$ )	Old BBN designation	IEST designation	ASHRAE designation	Typical Application
100 (4000)			Curve F	Surgical suites
50 (2000)	BBN-A	VC-A	Curve E	“Generic” lab
25 (1000)	BBN-B	VC-B	Curve D	Non-Photo Semiconductor
12.5 (500)	BBN-C	VC-C	Curve C	
6.3 (250)	BBN-D	VC-D	Curve B	Semiconductor Photolithography
3.1 (125)	BBN-E	VC-E	Curve A	Nanotechnology

Hal Amick & Michael Stead, “Vibration Sensitivity of a Laboratory Bench Microscope,” Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Variation of Sensitivity



Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Floor Criteria for Microscopes

Magnification	Benchttop	Newport Table
40x – 100x	100 $\mu\text{m/s}$ (4000 $\mu\text{in/s}$ )	1000 $\mu\text{m/s}$ (40,000 $\mu\text{in/s}$ )
400x	VC-A 50 $\mu\text{m/s}$ (2000 $\mu\text{in/s}$ )	500 $\mu\text{m/s}$ (20,000 $\mu\text{in/s}$ )
1000x	VC-C 12.5 $\mu\text{m/s}$ (500 $\mu\text{in/s}$ )	125 $\mu\text{m/s}$ (5000 $\mu\text{in/s}$ )

Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico

# Conclusions

- Microscope equally sensitive in all directions
- Greatest sensitivity between 10 and 50 Hz
- Sensitivity increases with magnification
  - Relatively insensitive at 40x
  - Order of magnitude increase in sensitivity between 40x and 1000x
- Optical bench has major effect (~20 dB)

Hal Amick & Michael Stead, "Vibration Sensitivity of a Laboratory Bench Microscope," Presented at the First Pan-American/Iberian Meeting on Acoustics; 144<sup>th</sup> Meeting of the Acoustical Society of America, 2-6 December 2002, Cancun, Mexico